

 CRISTALENS

CRISTAL<sup>®</sup>

HYDROPHILIC LENS



► Monofocal

MADE IN FRANCE

CE 0459

# TECHNICAL SPECIFICATIONS

<b>Lens type</b>	For implantation in the capsular bag
<b>Optic diameter</b>	6.00 mm
<b>Overall diameter</b>	13.00 mm
<b>Design</b>	One-piece square edge on 360°
<b>Optic design</b>	Monofocal Aspherical with negative spherical aberration to partly correct corneal spherical aberration Biconvex
<b>Angulation</b>	10°
<b>Material</b>	Hydrophilic Benz 25 made of hydrophilic and hydrophobic copolymers
<b>Dioptric powers</b>	From +10.0D to +30.0D by 0.5D
<b>Estimated A-Constant (SRK-T)</b>	118.0 Ultrasound biometry 118.5 Interference laser biometry
<b>Suggested anterior chamber depth (ACD)</b>	4.96 mm Ultrasound biometry 5.25 mm Interference laser biometry
<b>Refractive index</b>	1.46
<b>Sterilization</b>	Steam
<b>Recommended injection system</b>	Single-use injector
<b>Recommended incision size</b>	2.2 mm

