


ARTISAN[®] APHAKIA

#1 Backup Lens in complicated Cataract Cases


Based on the long term experience of Iris Fixation, the ARTISAN[®] Aphakia IOL is a predictable, safe, high precision implant, that corrects the eye when it is not correctable by other means. Pediatric lenses are available upon request for juvenile or congenital cataract cases.



Main Features & Benefits

- ✓ IRIS FIXATION
- ✓ IRIS CLAW[®] DESIGN
- ✓ LONG TERM CLINICAL EXPERIENCE
- ✓ PREDICTABLE , STABLE, RELIABLE

ARTISAN[®] Aphakia | Model 205



Optic:	5.0 mm
Haptics:	Iris Claw [®]
Overall Ø:	8.5 mm
A-Constant:	115.0
AC Depth:	3.3 mm
Dioptic Powers:	2.0 D to 30.0 D (1.0 increments) 14.5 D to 24.5 D (0.5 increments)
Also available:	ARTISAN [®] Pediatric Aphakia 4.4/6.5 and 4.4/7.5 Designed for small eyes

ARTISAN[®] Training Course

To obtain the best results, OPHTEC continuously conducts specialized training courses throughout the world to develop and maintain surgical skills. OPHTEC's condition to use the product correctly is to provide ARTISAN[®] IOLs only to certified professionals.

For more information on training courses send an e-mail to marketing@ophtec.com, or call +31 50 525 1944.